



Schure Foot Stainless Steel

Surgical positioning device for total knee replacement procedures

Easily mounts to operating room side rails using the Schure Socket XL, letting the surgeon replace the entire knee while the patient's foot is held firmly in place.

Features

- Disposable conductive, non-latex foot pad included
- Patient's well leg lies flat
- 12"L x 2"D, 2 lbs
- Schure Socket XL sold separately

#800-0062 Schure Foot, Stainless Steel

#508-0111 Disposable Foot Pads, 12/cs

Schure Socket XL

#800-0134 US

#800-0134-EU EU

#800-0134-JPN JPN

#800-0134-UK UK

#800-0134-DEN DEN

