



# Standard Armboard

*Offers a full range of motion, swings 180°*


Simply pull finger trigger at end of armboard and move to required position. Release trigger and it locks tight. To remove, just pull finger trigger and lift off side rail.

## Features

- Simply place armboard on OR table side rail and it snaps locked tight
- Pad sold separately
- 26.5"L x 6"W

## Standard Armboard

- |           |                       |           |     |
|-----------|-----------------------|-----------|-----|
| #800-0008 | US                    | #800-0189 | UK  |
| #800-0188 | EU                    | #800-0190 | DEN |
| #508-0071 | 2" (5cm) Deluxe Pad   |           |     |
| #508-0151 | 2" (5cm) Softcare Pad |           |     |



*New trigger for quick & easy removal*

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

452 Randolph Street, Abington, MA 02351

©2018 Schuerch Corporation #065-0047 REV 2.00