

Wrist Tower

Innovative Wrist Tower offers a wide variety of positioning capabilities for wrist procedures

Features

- *Simply dial in your tension...eliminates use of sand bags*
- *You can apply traction by simply dialing the knob from 0-25 lbs (0-11 kgs)*
- *Pronation and supination can be easily controlled by the surgeon's preference and locked in place*
- *Easily attaches to side rail of surgical table with CamLoc Clamp (sold separately)*
- *Improves productivity... eliminates the need for assistant to hold extremity*
- *All components are autoclavable*

#800-0135 Wrist Tower

Required Accessory: Order One CamLoc Clamp

#800-0271 US **#800-0271-UK** UK
#800-0271-EU EU **#800-0271-DEN** DEN

Optional Accessory

Finger Traps, Limited Reuse

Please call customer service for product pricing and specifications.

